

Introduction of IPK's Technology Platforms

2022 Pasteur Network
Annual Meeting


Regis Grailhe, PhD
& David Shum





Contents

- Drug screening & facilities
 - Protein-protein interaction studies
 - Technologies for resource limited environment
-



1 Target and drug screening facilities & technologies

David Shum,
Head of Screening
Discovery Platform
Seoul, South Korea



IP
IP Morocco
IP Montevideo

Biological Safety Level 2+ (Three Systems)

HCS: Two Confocal & Epi-Fluorescent >20,000 pts/day

HTS: >25,000 pts/day

BSL2

Biological Safety Level 3 (One System)

HCS: Epi-Fluorescent w/ Confocal mode >10,000 pts/day

HTS: >20,000 pts/day


BSL3



1 IPK Screening Technologies

Workflow & Access Points

Small Molecule & Natural Products & RNAi Libraries





Tools & Automated Technologies




Phenotypic & Target

iPSCs




Advanced Cellular Models & Assays




Spheroids & Organoids

Cell Painting & Omics



Ai-driven Analysis & Target ID



Model Prediction


Drug discovery pipeline




1

HTS technology: target-based

Biochemical model




Multimodal plate reader platforms



- 1 – known target essential for disease pathogenesis
- 2 – suitable for HTS
- 3 – capture single & multiple data point
- 4 – rapid analysis

1 HCS technology: image based cell-based Phenotypic-based

Cellular disease model



Phenotypic screening Automated microscopy platforms




- 1 – relevant models as close to physiological condition
- 2 – capture many types of parameters
- 3 – more information = High content


1

Image-based Profiling: Cell Painting

Parasite infected cell




Non-infected cell




Parasite infected cell


Hoechst33342
DNA (Nuclear stain)




Fluor488
ER (Concanavalin A)




Fluor512
RNA (Nucleic acid stain)




Fluor568
Actin (Phalloidin)



Fluor641
Mitochondria (MT stain)




Overlay




1 Activity Profiling: MOA & Target Deconvolution


1 Feature



~100 Features



Robust features (blue Low, Red High)




1 Libraries: starting points of discovery

Compound management

Drug: Chemical Libraries

- 200,000 small molecules
- 200,000 natural products




Target: high-throughput knockdown

- Arrayed siRNA library
- Arrayed lentivirus shRNA libraries




2 Screening technologies

Protein-protein interaction (BRET)




Study protein network

- Study the full combinatory of interaction of a subset of proteins at the cellular level.
- Determine the precise interaction location at sub-cellular level.




Fast kinetics (10Hz)

- Study the kinetic of interaction (ms, hour, or day scale).




Drug Screening

- Screening for PPI drug modulators.



Compare PPI relative affinity


- Determine the affinity of your protein against targets in living cells.
- Study effect of mutations.



2 Protein-protein interaction screening

Live cell PPI (BRET)

Dr **Regis Grailhe**,
Head of Technology
Development Platform
Seoul, South Korea




- *Live virus-host interaction characterization (Kinetics)*
- *Relative affinity (BRET50)*
- *Subcellular PPI imaging*

BRET nanoluciferase assay



Dr **Caroline Demeret**,
Head of Molecular
Genetics of RNA viruses
Paris, France



- *HTS capabilities*
- *Intensity map of interactions*
- *Comparative virus-host interaction profiling*

Split nanoluciferase assay



Acknowledgements

IPK Biology Teams

Translational Teams – SDP, MC (Inhee Choi) & TDP

Thank You

siRNA and drug screening

Contact: david.shum@ip-korea.org

Protein-protein interaction and low cost diagnostic

Contact: regis.Grailhe@ip-korea.org




For Research
For Health
For Our Future



3

Developing technologies for resource limited environment

Design and making diagnostic tools




Step 1: Contact regis.Grailhe@ip-korea.org

Step 2: Discussion (Skype, Microsoft Team, Zoom)

Step 3: Design

Step 4: Making and validation

Step 5: Data transferred by Email




3

Developing technologies for resource limited environment

Hardware development




Dr Regis Grailhe,
Head of Technology
Development Platform
Seoul, South Korea



Dr Lawrence Ayong,
Head for the Malaria
Research Laboratory
Yaounde, Cameroon



Imager*
bioluminescence



**Smartphone
microscope adapter**



LAMP (12V)
Loop mediated isothermal amplification

*doi.org/10.1021/acssensors.2c00457

